



**p/n 00-632 Suggested List: \$98.00 ! For limited time : With Travel Case**

# **PULSE OXIMETER**

This rugged and compact Finger-Tip Pulse Oximeter is a portable (includes lanyard), non-invasive blood oxygen saturation meter using two complimentary wavelengths of light for rapid and accurate %SpO<sub>2</sub> measurement. Display includes real-time Pulse Rate measurement with Pleth Waveform indication if desired.

Applicable for infants, children and adults, it can be used for reliable checks in settings as diverse as the hospital, office, clinic, aviation and field. Low power consumption, AAA batteries (included) and auto shut-off combine with a variable brightness Organic LED display for over 40 hours of continuous use.

- Certifications: FDA(USA), PMDA(Japan), CMDCAS(Canada), CE
- SpO<sub>2</sub>: 70-79% ( $\pm 3\%$ ), 80-99% ( $\pm 2\%$ ), Resolution  $\pm 1\%$
- Pulse Rate: 30-100 ( $\pm 2$ BPM), 101-235 ( $\pm 2\%$ ), Resolution  $\pm 1$  BPM
- Dual color OLED display, 10 Level Brightness
- Display SpO<sub>2</sub>, pulse rate, pulse bar and plethysmogram
- Six display modes, user selectable (see above images)
- Size: 58mm x 32mm x 34mm Weight : 28g (without 2xAAA batteries)
- Low power indicator (~40hrs), Auto-Off if no finger inserted
- One Year Warrantee (exclusive of abuse or damage)



Brochure 00-632.doc

## **Anesthesia Associates, Inc.**

460 Enterprise Street  
 San Marcos, CA 92078  
 Tel: (760) 744-6561 Fax: (760) 744-0054  
 Email: Solutions@AincA.com Web: www.AincA.com