

ADAPTERS AND CONNECTORS

CHROME-PLATED BRASS, (*LF)

Mask Adapter

#00-129

15mm male on the machine end. Patient end is nominal 22mm male and 15mm female.

#00-129



Hanks-Rackow Elbow

#00-132

15mm male on machine end. Patient end is nominal 22mm male and 15mm female.

#00-132



Hanks-Rackow Elbow with Side Inlet

#00-133

15mm male on machine end. Patient end is 22mm male and 15mm female with a 1/4" side inlet for introducing fresh gas at the patient.

#00-133



Universal Adapter

#00-320

Large end is 22mm female. Small end is 15mm female.

#00-320



Keets Type Hanks-Rackow Elbow, Side Inlet #00-443-A

15mm male on machine end. Patient end is 22mm male and 15mm female with a 1/4" side inlet which extends down to within 9/16" of the mask fitting, permitting an ideal "T" Tube principle for the flushing of the exhaled gas from under a mask.

#00-443-A



Keets Type Hands-Rackow Elbow, Top Inlet #00-443-B

15mm male on machine end. Patient end is 22mm male and 15mm female with a 1/4" top inlet which extends down to within 9/16" of the mask fitting, permitting an ideal "T" Tube principle for the flushing of the exhaled gas from under a mask.

#00-482

Keets Type Hanks-Rackow Elbow, Top Inlet, Side Port

#00-482

15mm male on machine end. Patient end is 22mm male and 15mm female with a 1/4" top inlet which extends down to within 9/16" of the mask fitting, permitting an ideal "T" Tube principle for the flushing of the exhaled gas from under a mask. Also has a 1/4" side port.



#00-443-B