

STETHOSCOPES AND ACCESSORIES

CHEST PIECES (*LF)

Wenger

Adult #00-390-A
Child #00-390-C

Heavyweight chest piece of chrome-plated brass. May be used on some areas of the body without being taped. Male Luer fitting for attachment.



Martin

#00-441

Lightweight Martin type, accommodates standard plastic I.V. tubing. The angle of the inlet keeps the chest piece flat. One size for both premature and infant. Chrome-plated brass.



Standard

Premature #00-442-P
Infant #00-442-I

Lightweight with male Luer fitting. Chrome-plated brass.

#00-442-I


#00-442-P



Standard, non-metallic

Premature #00-541
Child #00-542

Lightweight with male Luer fitting. Delrin/Nylon. MRI compatible.



Double-Sided Adhesive Disks

#00-390-DISC

Circular disks with center removed. Aggressive adhesive on both faces adhere securely any of the above chest pieces to the patient. Available in rolls of 100 qty, or by the 50 roll case (not shown).

MONOSCOPE (*LF)

		Left	Right
Monaural Stethoscope	Small	#00-383	#00-380
	Medium	#00-384	#00-381
	Large	#00-385	#00-382

Monaural stethoscope complete with PVC sound tubing and female Luer connector, comfortable single ear piece for constant monitoring of the patient.

The monaural stethoscope is molded of sufficient durometer to accurately reflect sound, yet soft enough to assure comfort while being worn. The ear piece is more or less C-shaped and is inserted into the grooves of the ear.

Clinical experience indicates that the medium size normally fits the greater proportion of users.

