

#00-204  
#00-218  
#00-220

# FLEX CONNECTORS COMPONENTS AND ADAPTERS

## TUBING (\*LF)

**Tubing (EPDM, reusable) (\*LF) #00-217**

EPDM synthetic rubber 5/8" tubing. Specify length up to 36".  
Uses adapter end fittings shown below.

## ADAPTERS (\*LF)

**Male Adapter** delrin #00-204  
15mm male x 5/8" flex stainless steel #00-218  
aluminum #00-220

**Female Adapter** delrin #00-205  
15mm female end x 5/8" flex stainless steel #00-219  
aluminum #00-221

**Tubing Adapter, Aluminum #00-481**  
22mm female end x 5/8" flex

## TUBING, EVACUATION (\*LF)

Flexible EPDM synthetic rubber tubing with integral 19mm female cuffs. Order by overall length as indicated.

**19mm female x 19mm female x 10" #LF-217-10**

**19mm female x 19mm female x 30" #LF-217-30**

**19mm female x 19mm female x 39" #LF-217-39**



#00-217



#00-205  
#00-219  
#00-221

#00-481



#LF-217

\*See page 43 for cleaning/sterilization information