



Our newly repackaged and upgraded Peripheral Nerve Stimulator is now available! It features:

- Fully adjustable output voltage
- Automatic feature repeat
- User friendly controls with indicator LEDs
- ALL functions are selectable from the front touch panel
- Sealed front membrane panel
- Use with ball electrodes or lead sets as required
- 9v battery with easy battery access
- Hand held, ergonomic design with non slip rubber sides
- Light weight yet rugged construction
- Proudly made in the USA!

Functions Include:

- Twitch of 1 per 8 seconds, repeating
- Train of Four (TOF), repeating
- 50Hz Tetanus
- 100HZ Tetanus
- ON and low battery LED indicator
- Pulse generated LED indicator
- Selected function LED indicators
- Adjustable output voltage control knob



Ordering Part Numbers

00-100-PLUS	Ainca Peripheral Nerve Stimulator PLUS with ball electrodes and 9v battery.
00-100-PLUS-KIT	The Ainca Peripheral Nerve Stimulator PLUS with ball electrodes, 9v battery, alligator clip lead set, and a soft sided case.

Accessories (ALL cross compatible with the 00-100)

00-106	Ball Electrode set
00-107-P	54" lead set with alligator clips
00-108-P	54" lead set with female snaps
00-109	Neoprene soft-sided zip case
00-110-P	54" extension lead set